


USE AN **AED** AS SOON AS POSSIBLE

AEDs are safe, easy to use, and can save a life.



1 **TURN ON**

the AED and follow the voice prompts.



2 **ATTACH PADS**

Place the pads on the person's bare chest as shown in the diagrams.



3 **CLEAR**

Make sure no one is touching the person when the AED analyzes.



4 **SHOCK (IF ADVISED)**

If the AED advises a shock, clear the person and press the shock button.

KEY DO'S & DON'TS

✓ **DO:**

- ✓ Push in the center of the chest.
- ✓ Push hard and fast.
- ✓ Allow full chest recoil.
- ✓ Switch with another person every ~2 minutes if possible.

✗ **DON'T:**

- ✗ Stop to check for a pulse.
 - ✗ Lean on the chest between pushes.
 - ✗ Give breaths if you're untrained.
- Hands-only CPR = compressions only.

WHEN TO USE HANDS-ONLY CPR



For adults and teens (ages 12 and older) who collapse suddenly.



Most cardiac arrests happen at home, at work, or in public places.



If you're untrained or unwilling to give rescue breaths, do hands-only CPR and call for help.



AED PURCHASE OPPORTUNITIES

Scan the QR code to find AED resources and purchase opportunities for your location.



911



REGISTER YOUR AED LOCATION

Register your AED with NH 911. Scan the QR code to register.



WATCH A SHORT HOW-TO CPR VIDEO

Scan the QR code to watch now.



CALL 911. PUSH HARD & FAST. USE AN AED.
YOU CAN SAVE A LIFE.